

Ultimative Speichelkontrolle

Weitaus beständiger als Watterollen aus Baumwolle und andere Saugeinlagen
Sie verlieren keine Flüssigkeit bei der Handhabung

Vorteile von DryTips® gegenüber herkömmlichen Watterollen aus Baumwolle:

- Effektive Abdeckung des Ohrspeicheldrüsen-Kanals und der Wangenschleimhaut
- Polyethylenfolie als Schnittschutz
- Erleichtert das Arbeiten weit hinten in der Mundhöhle
- Flexibel, passt sich Wangenbewegungen an
- Bleibt auf angenehme Weise am Platz und lässt sich gleichzeitig leicht mit Wasserspray entfernen
- Reflektierende Beschichtung für bessere Sicht im Arbeitsfeld



NEU!

Jetzt mit reflektierender Beschichtung erhältlich.

VERWENDUNG:


- Abformungen
- Okklusale Adjustierung
- Kieferorthopädie
- Fluoridanwendungen
- Komposite
- Zementierung
- Fissurenversiegelung
- Farbbeschichtungsverfahren

Für Ihre DryTips® Bestellung:

Artikelnummer	Beschreibung
291542	DryTips® groß, weiß, 50 Stk.
291543	DryTips® klein, weiß, 50 Stk.
291726	Reflective DryTips® groß, silber, 50 Stk.
291727	Reflective DryTips® klein, silber, 50 Stk.



Gebrauchsanweisung und Video unter www.microbrush.de



Ultimate moisture control

Far outlasts cotton rolls and other absorbent pads, and doesn't lose moisture when handled.

DryTips® advantage over traditional cotton rolls:

- Effective coverage of parotid orifice and buccal mucosa
- Polyethylen film for cutting protection
- Enables easier access to molar region
- Flexible, accomodates cheek movement
- Stays comfortably in place, yet is easily removed with water spray
- Reflective coating for increased visibility and to aid in retraction



NEW!
Reflective coating now available.

IDEAL FOR:

- Impressions
- Occlusal adjustments
- Orthodontics
- Fluoride applications
- Composites
- Cementation
- Fissure sealing
- Paint-on techniques

For your DryTips® order:

Item number	Description
291542	DryTips® large, white, 50 ct.
291543	DryTips® small, white, 50 ct.
291726	Reflective DryTips® large, silver, 50 ct.
291727	Reflective DryTips® small, silver, 50 ct.



Video and instructions for use available at www.microbrush.eu