


The Original Microbrush® Applicator

Consistent, reliable quality



NEW
DUO:
Double Bend

4

Tip sizes

ultrafine	0,5 mm
superfine	1.0 mm
fine	1.5 mm
regular	2.0 mm

11

Colours



3

Packaging options

Dispenser Kit
Dispenser Refill
Tube

Microbrush[®] DUO[™]




An Extension of You

Ergonomically designed, DUO's stay-put bend holds its shape so you can focus on the patient—moving seamlessly with you during every procedure.



Precision Redefined

DUO delivers exact, drip-free flow—applying the right amount every time for adhesives, etchants, sealants, and more, so you can work confidently without guessing.



Double Bend Design

DUO's dual bends reach tight spots others can't, and hold their shape, making precise application smooth and natural—giving you confidence and ease every time.

PRECISION, PERFECTED.

When your tools work as hard as you do, every procedure gets sharper, smoother, and more successful. Microbrush DUO redefines what's possible with a revolutionary double-bending tip that grants effortless access to the toughest areas, precision control that eliminates waste, and ergonomic comfort that feels like second nature. Engineered for professionals who demand consistency, performance, and reliability—DUO doesn't just keep up, it elevates. Ready to reach further, deliver smarter, and work more confidently?

DUO DELIVERS.



DOUBLE-BEND DESIGN

DUO's innovative double-bending tip unlocks access to the most challenging areas, giving clinicians precision where it matters most.



ULTIMATE CONTROL

Designed to feel like an extension of your hand, DUO offers unmatched stability and ergonomic ease during every procedure.



OPTIMAL DELIVERY

From etchants to adhesives, DUO ensures consistent, controlled application for flawless results every time.



Further information

www.microbrush.eu/en/microbrush-duo-en/

Item #	Color	Quantity
292724	Pearl Purple	100ct Tray
292725	Pearl Blue	100ct Tray

MB-EN-CAT/0126

We are happy to help!



info@microbrush.eu | + 49 (0) 6221 43 45 442

YOUNG[®]
INNOVATIONS

UltraBrush[®] Bristle Brush



- Bristle brush applicators for spreading and smoothing materials
- Bends close to tip for improved control
- “Stay-Put Bending” Technology - Maintain the applicator’s head at desired angle without rebound
- Chemical-resistant fibres designed to hold up to the harshest dental materials
- Available in two sizes: UltraBrush[®] regular ± 200 fibres, UltraBrush[®] fine ± 100 fibres
- Ideal for desensitizers, whitening solutions, fluoride, sealants and more



Item #	Description		
U2D	UltraBrush [®] Dispenser Kit	● green	regular
U2R200	UltraBrush [®] Refill	● yellow / green	regular
U1D	UltraBrush [®] Dispenser Kit	● blue	fine
U1R200	UltraBrush [®] Refill	● orange / blue	fine



Item #	Description			
MPD	Microbrush® Plus Dispenser Kit	● yellow	fine	1.5 mm
Microbrush® Plus Applicator Dispenser Kit				
PR400KIT	Microbrush® Plus Applicator Dispenser Kit	assorted	regular	2.0 mm
PF400KIT	Microbrush® Plus Applicator Dispenser Kit	● ● pink/yellow	fine	1.5 mm
PSF400KIT	Microbrush® Plus Applicator Dispenser Kit	● white	superfine	1.0 mm
Microbrush® Plus Applicator Refill				
PR400	Microbrush® Plus Applicator Refill	assorted	regular	2.0 mm
PR400BL	Microbrush® Plus Applicator Refill	● blue	regular	2.0 mm
PR400GR	Microbrush® Plus Applicator Refill	● green	regular	2.0 mm
PR400PE	Microbrush® Plus Applicator Refill	● peach	regular	2.0 mm
PR400PU	Microbrush® Plus Applicator Refill	● purple	regular	2.0 mm
PF400	Microbrush® Plus Applicator Refill	● ● pink/yellow	fine	1.5 mm
PF400PI	Microbrush® Plus Applicator Refill	● pink	fine	1.5 mm
PF400YE	Microbrush® Plus Applicator Refill	● yellow	fine	1.5 mm
PSF400	Microbrush® Plus Applicator Refill	● white	superfine	1.0 mm
PU400TE	Microbrush® Plus Applicator Refill	● teal	ultrafine	0.5 mm
PU400OR	Microbrush® Plus Applicator Refill	● orange	ultrafine	0.5 mm
Microbrush® Tube Applicators				
MRA400	Microbrush® Tube Applicators	assorted	regular	2.0 mm
MRB400	Microbrush® Tube Applicators	● blue	regular	2.0 mm
MRG400	Microbrush® Tube Applicators	● green	regular	2.0 mm
MRS400	Microbrush® Tube Applicators	● peach	regular	2.0 mm
MRP400	Microbrush® Tube Applicators	● purple	regular	2.0 mm
MFA400	Microbrush® Tube Applicators	● ● pink/yellow	fine	1.5 mm
MFP400	Microbrush® Tube Applicators	● pink	fine	1.5 mm
MFY400	Microbrush® Tube Applicators	● yellow	fine	1.5 mm
MSF400	Microbrush® Tube Applicators	● white	superfine	1.0 mm
MUT400	Microbrush® Tube Applicators	● teal	ultrafine	0.5 mm
MUO400	Microbrush® Tube Applicators	● orange	ultrafine	0.5 mm






Microbrush[®] X

- For even placement of material to inner surface of tooth wall
- Extended and extra thin for use in confined spaces
- Designed with stiffer fibres for scrubbing
- For endodontic treatment, also for fissures and bonding posts
- Patented bending portion is bendable up to 90 degrees



MPX
Microbrush[®] X Dispenser
Kit, black, 100 ct.

PX100
Microbrush[®] X Refill, black,
100 ct.





STIX64B
Micro-Stix™ Applicator,
64 ct., blue

Micro-Stix™ Applicator

- Ideal for veneers, crowns, brackets, inlays and onlays
- Applicator tip bendable
- Innovative adhesive formula
- Enables a secure hold of restorations until placed securely in their proper place.
- Can be used multiple times during a single patient procedure

